

In partnership with the Zomet Institute, AMIGO MOBILITY announces the...

# NEW Shabbat Amigo

- Shabbat "module" circuit board mounted inside throttle enclosure
  - Complex wiring eliminated
- ONE on/off key switch for either Shabbat/Normal mode

*Halachically-authorized for use on Shabbat by the Zomet Institute in Israel.  
Every Shabbat Amigo is individually inspected and certified by a Zomet representative.*

**NEW Shabbat module may be added to these Amigo models for an up-charge**  
*(in addition to the price of each Amigo.)*

- Amigo **RT Express** and **RT Express Junior** with 3" shorter platform option
- Amigo **RD** and **EXT** with 4" longer platform option
- Amigo **HD**

- 1 **On/Off Key Switch**  
(Normal or Shabbat mode)
- 2 **Shabbat/Normal Toggle Switch**
- 3 **Directional Rocker Switch**  
(Shabbat mode)
- 4 **Emergency Stop Rocker Switch**  
(Shabbat mode)
- 5 **Shabbat LED**



Amigo® RT Express



Amigo® RD



Amigo® HD



*Certification label appears on each Shabbat Amigo*

## SHABBAT CONTROL

Certified by

The **ZOMET Institute**

Approval #: \_\_\_\_\_

Type: \_\_\_\_\_

Model: \_\_\_\_\_

Serial #: \_\_\_\_\_



Your Amigo Representative is:

# SHABBAT SHALOM TO ALL!

Amigo Mobility International, Inc • 6693 Dixie Highway • Bridgeport, MI 48722-9725 • 1-800-821-2710 • www.myamigo.com  
Intended only for individuals whose limited mobility is dependent on a POV/scooter or automobile consistently throughout the week • M1386-09